


# ESSENTIAL DEVICES

## InSitu Accessory Portfolio

---

### Rainbow™ Angiography Catheters

- High level of torqueability
- Steel wire braided bodies
- Soft atraumatic tips
- Ergonomic hub
- Flexible tip segment
- High radiopacity
- Available for both Cardiology and Radiology shapes
- 4, 5, and 6F sizes
- 100 and 110 cm working lengths
- Multipacks available
- Judkins right, judkins left, pigtail, amplatz and more. For a complete list of tip varieties and ordering information contact sales



### StarFlex™ Sheath Introducer Kits

- Convenient and simple kit
- Two-way stop cock
- Arterial needle
- Low insertion force hemostasis valve
- Side port
- Mini guidewire
- Colored coded caps
- Individual sterile packs
- Pediatric 4 and 5F - 5.75 cm working lengths
- Femoral and transradial 4, 5, 6, 7, and 8F - 11/12 cm working length



### InSituSeal™ Hemostatic Patch System

- Hemostatic Gauze, becomes a viscous gel to accelerate micro-thrombus formation.
- Foam Backing, applied with the gauze to maintain a moist environment and allow easy pressure application.
- Foam Dressing, with adhesive, to maintain the hemostasis in a moist environment accelerating the healing.



## Intrac™ Endoflator

- 20cc reservoir
- 30 ATM
- Easy to use design
- Luminescent dial
- Ergonomic handle design
- Positive locking mechanism
- Quick latch for instant deflation
- Unobstructed view for bubble purging
- Angled gauge for easy reading
- Reinforced hose
- Clear barrel
- One handed release capabilities



## InSitu Guidewire Family

### Radial™ Interventional

- Excellent kink resistance
- Variety of tip flexibilities and configurations with Nitinol core
- Tested for multiple insertions
- Enhanced trackability with platinum distal tip
- Superior guidewire visualization
- Broad range of tip flexibility
- 0.014" diameter
- 190 and 300 cm lengths

### Radial™ Diagnostic Guidewires

- Enhanced trackability with teflon coated
- Superior guidewire visualization
- Available in fixed cores
- 0.021", 0.025", 0.035", and 0.038" diameters
- 150 and 260 cm working lengths



## Angioplasty Kit (PTCA Kit)

- Hemostasis Valve with extension and stop cock
- Torque device
- Insertion tool for PTCA wire



*The content in this brochure is provided for informational purposes only and is not intended or recommended as a substitute for educational professional medical advice or diagnosis purposes. The information is for use only in countries with applicable health authority product registration.*

INSITU TECHNOLOGIES® INC.  
539 Phalen Blvd, St. Paul, MN 55130, USA  
Tel: +1 651 389 1017  
Email: [sales@insitu-tech.com](mailto:sales@insitu-tech.com)  
Web: [www.insitu-tech.com](http://www.insitu-tech.com)

FDA Registered Facility  
EN ISO 13485:2016 certified  
For International Distribution. Rev.07/29/2021  
©2021 COPYRIGHT InSitu Technologies® Inc.

