


SIRO™ DrugEluting Stent


Coronary Sirolimus Drug Eluting Stent



- *Reliability and security unite*
- *Clinically proven results in patients worldwide**
- *SIRO™ combines our trusted Direct-Stent® platform with the clinically proven Sirolimus drug.*



*Clinical Evaluation Report - WR2-CER002-V600



From small vessels to long lesions, SIRO™ is the clinically proven D.E.S.

Technical Data

Specifications	Description
Stent Design	Hybrid Design - Open and closed cell
Stent Material	L605 Cobalt Chromium Alloy
Wall Thickness	0.08 mm / 0.003"
Strut Width	0.07 mm
Balloon System	Rx Rapid Exchange Balloon Catheter
Balloon Material	Low-Compliant – Nylon Blend
**Nominal Pressure	10 ATM
Rated Burst Pressure	16 ATM
*Minimum Guiding Catheter	5F (2.0 - 4.0mm Diameter); 6F (4.5 - 5.0mm Diameter)
Drug	Sirolimus Analogue
Dose	1.4 µg/mm ²
Polymer	Biodegradable Stable Polymer
Crimped Profile	0.9 - 1.2 mm
Tip Entry Profile	0.017"
Delivery Catheter Length	~ 142 cm

* See individual manufacturer for Guide Catheter Fr equivalent

** Assure full deployment of stent. Deployment pressures should be based on lesion characteristics.

Ordering Information

Stent Diameter	Stent Length (mm)	Ordering Number
2.25	08, 13, 18, 23, 28, 33, 38	<div style="display: flex; align-items: center; gap: 10px;"> <div style="border: 1px solid black; border-radius: 10px; padding: 5px 15px; font-weight: bold;">70</div> <div style="border: 1px solid black; border-radius: 10px; padding: 5px 15px; font-weight: bold;">Stent Diameter</div> <div style="border: 1px solid black; border-radius: 10px; padding: 5px 15px; font-weight: bold;">Stent Length</div> </div> <p>For example: Stent Diameter: 2.25 Stent Length: 08 => Ordering number: 7022508</p>
2.5	08, 13, 18, 23, 28, 33, 38, 43, 48	
2.75	08, 13, 18, 23, 28, 33, 38, 43, 48	
3.0	08, 13, 18, 23, 28, 33, 38, 43, 48	
3.5	08, 13, 18, 23, 28, 33, 38, 43, 48	
4.0	08, 13, 18, 23, 28, 33, 38	
4.5 **	13, 18, 23, 28, 33	
5.0**	13, 18, 23, 28, 33	

**Custom size, special order only.

The content in this brochure is provided for informational purposes only and is not intended or recommended as a substitute for educational professional medical advice or diagnosis purposes. The information is for use only in countries with applicable health authority product registration.

INSITU TECHNOLOGIES® INC.
 1301 Corporate Center Drive, Suite 160
 Eagan, MN, 55121, USA
 Text/WhatsApp +1 651-434-0062
 Email: hello@insitu-tech.com
 Web: www.insitu-tech.com

FDA Registered Facility
 ISO 13485:2016 certified
 For International Distribution
 ©2026 COPYRIGHT InSitu Technologies® Inc.