


CELOSIA

Covered Stent™

Celosia™ Coronary Covered Stent



- Mono layer stent graft with a U.S. patented intravascular hybrid open-closed cell design
- Clinically proven results in patients worldwide*
- ePTFE porous membrane

"Single-Center Non-Randomized Clinical Evaluation of the Safety and Reliability of the Celosia™ Covered Stent in Patients with Severe Coronary Occlusions" - Report on file



*Minimize Your Risk.
Let Us Cover You.*

Technical Data

Specifications	Description
Stent Design	Direct-Stent - Patented Hybrid Cell Design, Stent Ends Flarable
Stent Material	L605 Cobalt Chromium Alloy
Stent Covering Polymer	ePTFE Microporous Polymer, IND 90-120 µm
Wall Thickness	0.086 ±0.010 mm
Balloon System	Rx Rapid Exchange Balloon Catheter
Balloon Material	Semi-Compliant - Nylon / Pebax Blend
Optimal Deployment Pressure	10 ATM
**Nominal Pressure	10 ATM
Rated Burst Pressure	16 ATM (2.25-4.0 mm); 14 ATM (4.5-6.0 mm)
Shaft Size (Proximal/Distal)	2F / 2.7F
*Minimum Guiding Catheter	5F (2.25-3.0 mm); 6F (3.5-6.0 mm)
Crimped Profile	0.95-1.4 mm
Tip Entry Profile	0.016-0.017 inch
Recommended Guidewire	0.014 inch
Delivery Catheter Length	142 ±10 cm usable length
Stent Length	10-30 mm standard; 6.0 mm custom lengths per ordering table
Stent Ends Flaring	2.25, 2.5, 2.75 → 3.5 mm; 3.0-4.5 → 5 mm; 5.0, 6.0 → 7 mm
Radial Force	3-2.5 mm 1.3 N Force

**See individual manufacturer for Guide Catheter Fr equivalent.*

***Assure full deployment of stent. Deployment pressures should be based on lesion characteristics.*

Ordering Information

Stent Diameter (mm)	Stent Length (mm)
2.25	10, 13, 16, 19, 23, 26
2.5	10, 13, 16, 19, 23, 26, 30
2.75	10, 13, 16, 19, 23, 26, 30
3.0	10, 13, 16, 19, 23, 26, 30
3.5	10, 13, 16, 19, 23, 26, 30
4.0	10, 13, 16, 19, 23, 26, 30
4.5	15, 19, 23, 26, 30
5.0	17, 20, 27
6.0*	17, 20, 27, 30, 34, 37

Ordering Number

76 + Stent Diameter

+ Stent Length

Example: 2.25 × 10 mm

Ordering number: 7622510

*Custom Orders (6.0 mm diameter).

The content in this brochure is provided for informational purposes only and is not intended or recommended as a substitute for educational professional medical advice or diagnosis purposes. The information is for use only in countries with applicable health authority product registration.

1301 Corporate Center Drive, Suite 160
Eagan, 55121, Minnesota, United States of America
hello@insitu-tech.com
651-389-1017

FDA Registered Facility
ISO 13485:2016 certified
For International Distribution
©2026 COPYRIGHT InSitu Technologies Inc.

